

Development of Iron Thin Films for Polarization Analysis of Ultracold Neutrons

2021-1-12

Hiroaki Akatsuka, Masahiro Hino, TUCAN collaboration

Contents

Introduction

- Purpose of the project
- Principle of UCN polarization measurement

Development and Evaluation of UCN polarization analyzer films

- Development of UCN polarization analyzer films
- B-H curve measurement by Vibrating Sample Magnetometry (VSM)
- Neutron reflectivity measurement

Conclusions and Outlook

Purpose of the project

Context

- UCN polarization films were tested in 2018, but did not achieve a good performance (polarization: 60% cf. Sean Hansen-Romu 2019)

Purpose

- Establish in-house production procedures of iron films for reliable preparation of high-performance films (in collaboration with Prof. M. Hino (KURNS, Kyoto Univ.))
- In addition to conventional aluminum substrate, investigate silicon substrate iron films
- Characterize the samples with complementary methods:
Vibrating sample magnetometry, cold neutron reflectometry (for Si substrate), (-> this talk), UCN (-> next talk)

	Advantages	Disadvantages
Al substrate	<ul style="list-style-type: none">• Results used in previous research• Can be made thin (~ 20 um) (small absorption)	<ul style="list-style-type: none">• Larger magnetic field required for saturation• Neutron reflectivity measurement not possible
Si substrate	<ul style="list-style-type: none">• Neutron reflectivity measurement possible• Smaller magnetic field required for saturation	<ul style="list-style-type: none">• Cannot be made thinner than < 0.1 mm (Absorption is about~ 2% larger than aluminum)


Principle of UCN spin analyzer

- UCN are polarized by a magnetized iron film
- The effective potential experienced by UCN:

$$V = V_F \pm \mu_n B = 209 \text{ neV} \mp 60.3 \text{ neV/T} \cdot B$$

- Fermi potential of Fe: $V_F \sim 209 \text{ neV}$
- UCN kinetic energy $\lesssim 300 \text{ neV}$

- **With $\sim 2 \text{ T}$ magnetization**
 - full separation of the UCN spin states
 - the film can be used **as a spin analyzer**



Simultaneous Spin Analyzer (SSA)


- **Simultaneous Spin Analyzer (SSA)**

- Polarization analyzer
- Thin Fe film in a permanent magnet: allows only a specific spin state to transmit
- Simultaneous measurement of UCN in each spin state
 - Selection of the spin state by RF spin flipper
- Measure the polarization from the number of UCN for each spin.

- **Requirements on the polarization films for the SSA**

- Large saturation magnetization ($\sim 2\text{ T}$)
(for high efficiency of spin analysis)
- Saturate with a low magnetic field ($\lesssim 10\text{ mT}$)
 - Low leakage field
 - Compact device
- Small absorption of UCN

S. Afach, et al, Euro. Phys. Jour. A 51, 143 (2015)



Development of polarization analysis film

Iron thin films were prepared using an ion beam sputtering system (IBS) at KURNS






Produced sample

Iron thin film

- thickness : 30, 50, 90 nm

Substrates


- Si (thickness $380\ \mu\text{m}$)
- Al (thickness $25\ \mu\text{m}$)



Substrate holder



Fe target

Ni/C target
(not used this time)



M. Hino et al, Nucl. Inst. Meth. A 797, 265 (2015)

BH curve measurement



H_c : Intersection of BH curve and $B=0$, H_s :H corresponding to 95% of B_{sat}

sample	H_c (Oe)	H_s (Oe)	$\mu_0 H_s$ (mT)
Si+Fe(30 nm)	16.5	29.8	2.98
Si+Fe(50 nm)	16.7	24.5	2.45
Si+Fe(90 nm)	20.3	45.5	4.55


sample	H_c (Oe)	H_s (Oe)	$\mu_0 H_s$ (mT)
Al+Fe(30 nm)	45.0	159	15.9
Al+Fe(90 nm)	25.8	134	13.4

- BH curves of each sample were measured by Vibrating Sample Magnetometry (VSM)
- **Samples saturates with low magnetic fields (requirements : $\mu_0 H_s \lesssim 10 \text{ mT}$)**
 - $\mu_0 H_s \lesssim 5 \text{ mT}$ (Si substrate)
 - $\mu_0 H_s \lesssim 15 \text{ mT}$ (Al substarte)

Comparison to previous TUCAN films

Sample	H _c (Oe)	H _s (Oe)
Al+Fe(30 nm) KURRI	25.8	29.8
Al+Fe(90 nm) KURRI	20.30	45.5
Al+Fe(400 nm) LeBow	33	161
Al+Fe(150 nm) Movatec	113	219

H_c: Intersection of BH curve and B=0,
H_s:H corresponding to 95% of B_{sat}



KURNS films saturated by a comparable field with LeBow films and much smaller magnetic field than the Movatec films.


Thanks to Blair and Beatrice for sending the LeBow and Movatec foils!

Comments (Hino-san):

- Higher H_c of 30 nm than 90 nm is within lot-per-lot fluctuations. This level of fluctuations exist among different lots, but the quality is rather homogeneous in the same lot (see p.7)

Comparison to PSI films


- PSI film VSM measurement:
[annual report 2012](#)
- Al-substrate films by KURNS
are comparable to the PSI's.




Neutron reflectometry

- Cold neutrons (wavelength $0.2 \sim 1$ nm) are injected to a sample and the reflectivity is measured.
- From the reflectivity profile as a function of wave vector transfer $q = 4\pi \sin \theta_0 / \lambda$, the critical value q_c can be determined.
- From the critical value q_c , the magnetic potential experienced by the neutron can be extracted.


$$V_{\text{eff}} = \frac{\hbar^2 q^2}{8m_n} \quad (m_n : \text{neutron mass})$$



Journal of the Neutron Science Society of Japan "Ripples" Vol.18, No.4, 2008 Principle of neutron reflectometry Naoya Torikai and Masayasu Taketa



Neutron reflectivity measurement Setup (J-PARC/MLF BL05)








Reflectivity measurement of the sample

Measurement procedure:

1. A magnetic field of -8 mT was applied to the sample. (**Reset**)
2. The reflectivity R was measured while increasing the applied magnetic field.
 - The effective potential was determined from the reflectivity with the spin flipper on and off.

How to extract physical properties of interest:

- The critical value q_c of the reflectivity profile $R(q)$
 → The potential V_{eff} experienced by neutrons at saturation
- From the magnetic field at which the effective potential V_{eff} rises
 → Magnetic field required to saturate the film



Neutron reflectivity measurement results

Fe 90 nm + Si substrate
(Neutrons incident in each spin states
Reflectivity R_{down} , R_{up})

- Fitted with a function which takes into account the beam polarization
- Fitting results

$$V_{\text{eff}} = 308(1) \text{ neV}$$


(requirements: $V_{\text{eff}} > 300 \text{ neV} \gtrsim E_{\text{UCN}}$)

\leftrightarrow Fe magnetization 2.02(1) T (E_{UCN} : UCN energy)

when 8 mT is applied

(requirements : $\mu_0 H \lesssim 10 \text{ mT}$)

- **Enough saturation magnetization obtained to polarize UCN at sufficiently low magnetic field**



wave vector transfer q vs
Spin-up reflectivity R_{up} , Spin-down reflectivity R_{down}
when 8mT is applied.

The global fit was performed using the Fermi potential of iron and the magnitude of the saturation magnetization as parameters.

Conclusion and outlook

- Development of UCN polarization analyzer for TUCAN
- Reflectivity measurement of the sample (Fe film and Si substrate)
 - $V_{\text{eff}} \sim 308(1) \text{ neV}$ (requirements: $V_{\text{eff}} \gtrsim 300 \text{ neV} \gtrsim E_{\text{UCN}}$),
when 8 mT is applied. (requirements: $\mu_0 H \lesssim 10 \text{ mT}$)
- **Successfully produced films that can be used for UCN polarization analysis**
- **The iron thin film on Si substrate was found to be magnetized at about 1/3 of the applied magnetic field compared to the PSI thin film**
- **Evaluation with UCNs are planned in spring 2022.**

Acknowledgements

Thank you for your attention!



TUCAN
TRIUMF Ultra Cold
Advanced
Neutron source

